



Thank you for purchasing this Asanté IntraCore switch. This release note applies only to:

A. IntraCore 35160-T

Part number: 99-00747-xx
Sub-assembly number: 61-20747-00

B. IntraCore 35160-G

Part number: 99-00716-xx
Sub-assembly number: 61-20716-00

Version History

- 1.0 November 2003: Original release (CD Version 1.0)
- 2.01 October 2003: Bug fixes. Obsoletes 1.2
- 2.01 February 2004: Corrected username information

Factory Defaults

- *Console Port, SSH or Telnet*
User name: root
Password: Asante

- *Web Interface*
IP address is: 192.168.0.1
User name: IntraCore
Password: Asante

All entries are case-sensitive. It's **Asante** not asante.

Bug Fixes

1. If IGMP query is enabled, system will hang with machine check during the boot-up process.
2. IGMP query packet is changed to V1 general query type
3. Cannot create more than 1 trunk via web.
4. QoS display

Technical Support

For technical assistance or to download the latest firmware, visit www.asante.com/support